

NANOPARTZ™

The gold nanoparticle for nanotechnology

Certificate of Analysis

Product Family	Nanopartz Gold Nanorods			
Product Number	A12-10-1200-CTAB-DIH-1-5			
Lot #	12G126-12948			
Form	Gold Colloid			
Capping Agent	CTAB			
Solution	18MEG DI Water			
Source	12G126			
SPR Peak	nm	measured	UV/VIS	1200
SPR Abs	OD	measured	UV/VIS	1
LSPR Peak	nm	measured	UV/VIS	504
LSPR Abs	OD	measured	UV/VIS	0.25
pH		measured	Expt.	7
Zeta Potential	mV	measured	DLS	40
Diameter	nm	measured	TEM	10
Length	nm	measured	TEM	81
Aspect Ratio		calc		8.1
Concentration	nps/mL	calc		3.0E+11
Wt. conc.	mg/mL	calc		0.035
Concentration	ppm	calc		35
Wt %	%	calc		0.0035%
Molarity	nM	calc		0.496
Molar Ext.	M-1cm-1	calc		2.0E+09
CTAB	mM	measured		3
Volume	mL	measured		5
Date				7/8/2023
Exp Date				1/4/2024

Alex Schoen, Manager Quality Control

