

NANOPARTZ™

The gold nanoparticle for nanotechnology

Certificate of Analysis

Product Family	Accurate Spherical Gold Nanoparticles
Product Number	A11-15-CIT-DIH-1-25
Lot #	11N580-13603
Form	Gold Colloid
Capping Agent	Citrate
Solution	18MEG DI Water
Source	11H80

Diameter	nm	measured average	TEM, UVVIS	15
SPR Abs	OD	measured	UV/VIS	1
PDI	stddev/size	measured	TEM	10%
SPR Peak	nm	measured	UV/VIS	518
pH		measured	Orion	7
Zeta Potential	mV	measured	DLS	-33
Concentration	nps/mL	calc		1.47E+12
Wt. conc.	mg/mL	calc		0.05
ppm		calc		50
Wt. %	%	calc		0.0050%
Molarity	uM	calc		2.45E-03
Molar Ext.	M-1cm-1	calc		4.08E+08
Citrate	mM	measured		3
Volume	mL	measured		25

Date	1/11/2024
Exp Date	7/9/2024

Alex Schoen, Manager
Quality Control

