

# NANOPARTZ™

The gold nanoparticle for nanotechnology

## Certificate of Analysis

Product Family	Accurate Spherical Gold Nanoparticles
Product Number	A11-25-CIT-DIH-1-25
Lot #	11M551-13161
Form	Gold Colloid
Capping Agent	Citrate
Solution	18MEG DI Water
Source	11K329

Diameter	nm	measured average	TEM, UVVIS	25
SPR Abs	OD	measured	UV/VIS	1.1
PDI	stddev/size	measured	TEM	7%
SPR Peak	nm	measured	UV/VIS	521
pH		measured	Orion	7
Zeta Potential	mV	measured	DLS	-30
Concentration	nps/mL	calc		3.49E+11
Wt. conc.	mg/mL	calc		0.06
ppm		calc		55
Wt. %	%	calc		0.0055%
Molarity	uM	calc		5.82E-04
Molar Ext.	M-1cm-1	calc		1.89E+09
Citrate	mM	measured		3
Volume	mL	measured		25

Date	1/11/2024
Exp Date	7/9/2024

Alex Schoen, Manager  
Quality Control

