

NANOPARTZ™

The gold nanoparticle for nanotechnology

Certificate of Analysis

Product Family	Accurate Spherical Gold Nanoparticles
Product Number	A11-30-CIT-DIH-1-25
Lot #	11L494-13750
Form	Gold Colloid
Capping Agent	Citrate
Solution	18MEG DI water
Source	11J266

Diameter	nm	measured average	TEM, DLS, UVVIS	30
SPR Abs	OD	measured	UV/VIS	1.1
PDI	stddev/size	measured	TEM	9%
SPR Peak	nm	measured	UV/VIS	523
pH		measured	Orion	7
Zeta Potential	mV	measured	DLS	-31
Concentration	nps/mL	calc		2.02E+11
Wt. conc.	mg/mL	calc		0.06
ppm		calc		55
Wt. %	%	calc		0.0055%
Molarity	uM	calc		3.37E-04
Molar Ext.	M-1cm-1	calc		3.26E+09
Citrate	mM	measured		5
Volume	mL	measured		25

Date	1/11/2024
Exp Date	7/9/2024

Alex Schoen, Manager Quality Control

