

# NANOPARTZ™

The gold nanoparticle for nanotechnology

## Certificate of Analysis

Product Family	Accurate Spherical Gold Nanoparticles
Product Number	A11-75-CIT-DIH-1-25
Lot #	11H137-13398
Form	Gold Colloid
Capping Agent	Citrate
Solution	18 MEG DI water
Source	11H137

Diameter	nm	measured average	TEM, DLS, UVVIS	75
SPR Abs	OD	measured	UV/VIS	1.1
PDI	stddev/size	measured	TEM	4%
SPR Peak	nm	measured	UV/VIS	543
pH		measured	Orion	7
Zeta Potential	mV	measured	DLS	-35
Concentration	nps/mL	calc		1.29E+10
Wt. conc.	mg/mL	calc		0.06
ppm		calc		55
Wt. %	%	calc		0.0055%
Molarity	uM	calc		2.16E-05
Molar Ext.	M-1cm-1	calc		5.10E+10
Citrate	mM	measured		3
Volume	mL	measured		25

Date	1/11/2024
Exp Date	7/9/2024

Alex Schoen, Manager Quality Control

Absorbance

