

NANOPARTz™

The gold nanoparticle for nanotechnology

Certificate of Analysis

Product Family	Nanopartz Gold Nanorod
Product Number	A12-10-808-CTAB-DIH-100
Lot #	12G97
Form	Gold Colloid
Capping Agent	CTAB
Solution	18MEG DI Water

SPR Peak	nm	measured	UV/VIS	800
SPR Abs	OD	measured	UV/VIS	1
LSPR Peak	nm	measured	UV/VIS	511
LSPR Abs	OD	measured	UV/VIS	0.3
pH		measured	Expt.	6
Zeta Potential	mV	measured	DLS	35
Diameter	nm	measured	TEM	10
Length	nm	measured	TEM	40
Aspect Ratio		calc		4.0
Concentration	nps/mL	calc		6.3E+11
Concentration	mg/mL	calc		0.035
Concentration	ppm	calc		35
CTAB	mM	measured		3
Volume	mL	measured		100

Date	1/4/2017
Exp Date	7/3/2017

Alex Schoen, Manager Quality Control

