

NANOPARTz™

The gold nanoparticle for nanotechnology

Certificate of Analysis

Product Family	Accurate Spherical Gold Nanoparticles
Product Number	A11-40-CIT-DIH-1-10
Lot #	SPEC
Form	Gold Colloid
Capping Agent	Citrate
Solution	18MEG DI water
Source	11H104

Diameter	nm	measured average	TEM, DLS, UVVIS	40
SPR Abs	OD	measured	UV/VIS	1.1
PDI	stddev/size	measured	TEM	5%
SPR Peak	nm	measured	UV/VIS	527
pH		measured	Orion	7
Zeta Potential	mV	measured	DLS	-36
Concentration	nps/mL	calc		8.53E+10
Wt. conc.	mg/mL	calc		0.06
ppm		calc		55
Wt. %	%	calc		0.0055%
Molarity	uM	calc		1.42E-04
Molar Ext.	M-1cm-1	calc		7.74E+09
Citrate	mM	measured		5
Volume	mL	measured		10

Date 10/4/2019
Exp Date 4/1/2020

Alex Schoen, Manager Quality Control

