

# NANOPARTz™

The gold nanoparticle for nanotechnology

## Certificate of Analysis

Product Family	Accurate Spherical Gold Nanoparticles
Product Number	A11-85-CIT-DIH-1-10
Lot #	SPEC
Form	Gold Colloid
Capping Agent	Nanopartz Carboxylic Acid
Solution	18 MEG DI water
Source	11H137

Diameter	nm	measured average	TEM, DLS, UVVIS	85
SPR Abs	OD	measured	UV/VIS	1.1
PDI	stddev/size	measured	TEM	4%
SPR Peak	nm	measured	UV/VIS	558
pH		measured	Orion	7
Zeta Potential	mV	measured	DLS	-35
Concentration	nps/mL	calc		8.89E+09
Wt. conc.	mg/mL	calc		0.06
ppm		calc		55
Wt. %	%	calc		0.0055%
Molarity	uM	calc		1.48E-05
Molar Ext.	M-1cm-1	calc		7.42E+10
Citrate	mM	measured		5
Volume	mL	measure		10

Date	10/4/2019
Exp Date	4/1/2020

Alex Schoen, Manager Quality Control

