

NANOPARTz™

The gold nanoparticle for nanotechnology

Certificate of Analysis

Product Family	Nanopartz Gold Nanowires			
Product Number	A14-10000-CTAB-DIH-1-5			
Lot #	F3280-7999			
Appearance	Slight off white color			
Form	White precipitate colloid			
Capping Agent	CTAB			
Solution	18MEG DI Water			
Source	F3280			

SPR Peak	nm	measured	UV/VIS	510
SPR Abs	OD	measured	UV/VIS	0.05
pH		measured	Expt.	7
Zeta Potential	mV	measured	DLS	35
Diameter	nm	measured	TEM	75
Length	nm	measured	TEM	10000
Aspect Ratio		calc		133.3
Concentration	nps/mL	calc		5.9E+07
Wt. conc.	mg/mL	calc		0.050
Concentration	ppm	calc		50
Wt %	%	calc		0.0050%
Molarity	nM	calc		9.81E-05
Molar Ext.	M-1cm-1	calc		5.1E+11
CTAB	mM	measured		25
Volume	mL	measured		5

Date	8/28/2019
Exp Date	2/24/2020

Alex Schoen, Manager
Quality Control

