

Name	Abbreviation	Structure	Polarity	Acetic Acid	Acetone	acetonitrile	benzene	n-butanol	butyl acetate	carbon tetrachloride	chloroform	cyclohexane	1,2-dichloroethane	dichloromethane	dimethylformamide	dimethyl sulfoxide	dioxane	ethanol	ethyl acetate	diethyl ether	heptane	hexane	methanol	methyl-t-butyl ether	2-butanone	pentane	n-propanol	isopropanol	diisopropylether	tetrahydrofuran	toluene	trichlorethylene	water	xylene		
Acetic Acid	AcOH	C2H4O2																																		
Acetone	Ace	(CH3)2CO																																		
acetonitrile	ACN	CH3CN																																		
benzene	BEN	C6H6	Non																																	
n-butanol	n-BuOH	C4H9OH	Non																																	
butyl acetate	BUAC	C6H12O2	Non																																	
carbon tetrachloride	CTC	CCl4	Non																																	
chloroform	CHL	CHCl3	Non																																	
cyclohexane	Cy	C6H12	Non																																	
ethylene dichloride	EDC	C2H4Cl2	Non																																	
dichloromethane	DCM	CH2Cl2	Non																																	
dimethylformamide	DMF	(CH3)2NC(O)H																																		
dimethyl sulfoxide	DMSO	(CH3)2SO																																		
dioxane	DI	C4H8O2																																		
ethanol	EtOH	C2H5OH																																		
ethyl acetate	ETAC	CH3COOC2H5																																		
diethyl ether	Et	(C2H5)2O	Non																																	
heptane	NC7	C7H16	Non																																	
hexane	HX	C6H14	Non																																	
methanol	MeOH	CH3OH																																		
methyl-t-butyl ether	MTBE	C5H12O	Non																																	
Ethylmethylketone	MEK	CH3C(O)CH2CH3	Non																																	
pentane	IUPAC	C5H12	Non																																	
n-propanol	PrOH	CH3CH2CH2OH																																		
isopropanol	IPA	C3H8O																																		
diisopropylether	DIPE	C6H14O	Non																																	
tetrahydrofuran	THF	(CH2)4O																																		
toluene	MePh	C7H8	Non																																	
trichlorethylene	TCE	C2HCl3	Non																																	
water	DIH	H2O																																		
xylene	X	(CH3)2C6H4	Non																																	